

Aim :-


To identify the extra element and the functional group present in the given sample (salicylic acid).

Reference :-Requirement :-

Ignition tube, Pair of tongs, test tube holder, China dish, tripod stand, wire gauze, funnel and filter paper.

Theory :-Sodium bicarbonate test :


Place a little of the substance in a test tube and add 2 ml of 5% aqueous sodium bicarbonate solution. The acids form their corresponding sodium salts with the liberation of carbon dioxide.



Teacher's Signature _____

Ester test:-

Heat about 0.5 g of the sample with 1 ml of methyl alt alcohol and a few drops of concentrated H_2SO_4 in a dry test tube for about a minute. Cool and pour the contents to a few ml of water contained in a small beaker. Salicylic acid can be identified by the colour of wintergreen oil.



Concentrated H_2SO_4 checks the reversibility of the reaction and acts as a dehydrating agent.

Acriflavine test:

Appearance of yellow or brown yellow precipitate confirms the presence of $COOH$ group.

Ferric chloride test:

This test is carried out in a ~~control~~ neutral solution. Put about 0.5 g of the substance in a boiling test tube and add dilute ammonia solution until the contents are just alkaline.


Tests for Carboxylic Acids

S.No	Experiment	Observation	Inference
1.	Sodium bicarbonate test - Place a little of the substance in a test tube and add 2 ml of 5% aqueous sodium bicarbonate sol ⁿ .		
2.	Ester formation - Heat about 0.5g of the sample with 1 ml of methyl alcohol and add a few drops of conc. H ₂ SO ₄ in a dry test tube for about a minute. Cool and pour the contents to a few ml of water contained in a small beaker.		
3.	Acriflavine test. - To a mixture of 0.1% acriflavine and 1% potassium chromate solution add a few drops of aqueous solution of the given sample.		
4.	Ferric Chloride Test - Put about 0.5g of the substance in a boiling test tube and add dilute ammonia solution until the contents are just alkaline to litmus paper. Boil the solution gently until the odour of ammonia is not there. To cold solution add few drops of ferric chloride solution.		
5.	Benzylthiousonium test - Dissolve about 0.5g of the sample in minimum volume of water. Add 5% sodium hydroxide sol ⁿ until the solution is just alkaline to methyl orange. Add a drop of dilute HCl. In another test tube dissolve 0.8g of benzylthiousonium chloride in 2 ml of water. Mix the 2 solutions. Cool in ice bath.		

to litmus paper. Boil the solution gently until the odour of ammonia is ~~not~~ not there. To the cold solution add few drops of ferric chloride solution. A violet colour indicates the presence of salicylic acid.

Benzylthiuronium Test

Benzylthiuronium salt of the acid ~~separates~~ separates out. (Chemistry same as previous experiment).



Salicylic
acid

Ferric
chloride

Violet coloured
complex

Hydrochloric
acid.

Result:-

The given organic compound has been studied successfully in the laboratory.